PAPER • OPEN ACCESS


Infrared spectroscopy for a quantitative determination of $CH_{3-n}CI_nCOCI$ in $TiCI_4$

To cite this article: B Yang et al 2015 IOP Conf. Ser.: Mater. Sci. Eng. 103 012027

View the article online for updates and enhancements.

You may also like

- Ab initio calculations of nuclear guadrupole resonance frequencies in trichloroacetyl halides: a comparison of DFT and experimental data F Iwase
- <u>Isotopically selective dissociation of COCl</u> <u>molecules in the radiation field of a pulsed</u> <u>CO laser</u> T G Abzianidze, V Yu Baranov, A B Bakhtadze et al.
- <u>Synthesis, characterization and</u> <u>electrochemical properties of 5-</u> <u>aza[5]helicene-CH_O-CO-MWCNTs</u> <u>nanocomposite</u> F Fontana, F Melone, D lannazzo et al.





DISCOVER how sustainability intersects with electrochemistry & solid state science research



This content was downloaded from IP address 52.14.249.197 on 22/05/2024 at 00:14

Infrared spectroscopy for a quantitative determination of CH_{3-n}Cl_nCOCl in TiCl₄

B Yang, Y Wang, H Tan¹, J He and G Song

Guizhou Academy of Testing and Analysis, Guiyang 550002, China

E-mail: 357688237@qq.com

Abstract. Infrared spectroscopy permits a determination of the concentration of CHCl₂COCl, CH₂ClCOCl and CCl₃COCl as the organic impurities in TiCl₄In this study, an infrared spectrometer with ZnSe window and the cell of PTFE was assembled, and used to determine the concentration of three organic impurities. The detection limits of CCl₃COCl, CHCl₂COCl, CH₂ClCOCl were determined as 3.16×10^{-3} µg/g 1.92×10^{-3} µg/g and 1.55×10^{-2} µg/g respectively. The new method is applied to quality control of refined TiCl₄ from titanium sponge.

1. Introduction

Titanium is a transition metal with various excellent properties, such as low density and high corrosion resistance. Titanium sponge is important for titanium manufacture [1]. An enrichment correlation presented among the impurities in Titanium sponge and that in fine TiCl₄. Therefore, the impurities of TiCl₄ are the key of titanium quality control [2]. The composition and content of impurities are analyzed quickly and quantitatively to removal of impurities CH_{3-n}Cl_nCOCl (CCl₃COCl, CHCl₂COCl, CH₂ClCOCl) during the production of Titanium sponge. The quality control standard of TiCl₄ is strict in Japan, America and Russia, organic compound was determined by infrared spectral analysis, besides the determination of inorganic compound like As, Sb, Si, Cl, Al, Cu, Fe, Pb, V, Ni and Sn [3-6]. At present, the general determination of impurities in domestic titanium industry just includes SiCl₄, FeCl₃, VOCl₃, and chromaticity analysis [7]. The lack of effective determination method of organic impurities (CCl₃COCl, CHCl₂COCl, CH₂ClCOCl) reduces the products' hardness. However, CCl₃COCl in fine TiCl₄ was determined by infrared semi-quantitative method [8-9] and Fouriertransform infrared spectroscopy (FT-IR) [10].

SiCl₄, Fe Cl₃ and VOCl₃ are the main impurities in fine TiCl₄. It did not have any other interference to the determination of CH_{3-n}Cl_nCOCl due to they are inorganic compounds. FT-IR cannot detect the peak of Fe and Si since the different bond vibration compared to organic compound. Since TiCl₄ and CH_{3-n}Cl_nCOCl undergone rapid hydrolysis in the presence of moisture, an infrared liquid cell is damaged, and the results can be influenced. Therefore, a new liquid cell is designed with good tightness and adaptive thickness. Many other IR transparent inert solid window materials for liquid cells are known – some are cheaper and established. A superior waterproofness is the reason for the choice of ZnSe. It is safe, simple, easy removable, washable and reusable. Spiking method or standard addition [11] was done to determine the trace organic impurities due to the content of organic

Content from this work may be used under the terms of the Creative Commons Attribution 3.0 licence. Any further distribution (cc) $(\mathbf{\hat{H}})$ of this work must maintain attribution to the author(s) and the title of the work, journal citation and DOI. Published under licence by IOP Publishing Ltd 1

¹ Address for correspondence: H Tan, Guizhou Academy of Testing and Analysis, Guiyang 550002, China. Email: 357688237@qq.com.

impurities are quite low in TiCl₄.

FT-IR quantitative analysis was based on absorption strength and it was able to quantify the gas, solid and liquid sample. In the same time, the component of mixture was calculated through the characteristic peak. Therefore, it was able to determine the infrared spectrum of the mixtures by FT-IR without pre-processing. And this measurement method was convenient, visible and quickly with several calculation method. In the study, the component of fine $TiCl_4$ was determined by the Standard addition method.

2. Experimental procedure

2.1. Instrument and reagents

Reagents Titanium tetrachloride was purchased from Zunyi Titanium Company (Guizhou, China), Chloroacetyl chloride, Dichloroacetyl chloride and trichloro-acetic chloride were purchased from Sigma Chemical Co. (St. Louis, MO, USA). Aladdin, dry calcium chloride, absolute ethyl alcohol are purchased from Shanghai Chemical Company.

A Thermo (IR 200) FT-IR was used in the experiment, and it connected with a spectrophotometer operated at a resolution of 4 cm⁻¹. Thirty-two scans were scanning to improve the signal-to-noise ratio.

2.2. IR liquid cell

Make sealed IR liquid cell for avoiding TiCl₄ reacts vigorously with moisture. It is made of a Teflon body, a window of ZnSe (25 mm \times 2 mm, 7800 \sim 440 cm⁻¹) and soft Teflon gaskets. Filling the TiCl₄ from these cells into the atmosphere was also found to be a troublesome process for routine analysis. Hence, this IR cell was fabricated. Structure and practicality chart of IR seal cell showed below (figure 1).




Figure 1. Structure and practicality chart of assembled infrared seal pool.

2.3. Sample preparation

Organic impurities mixture standard stock: Addition of 5 μ L CCl₃COCl (98%), 9 μ L CHCl₂COCL (98%), 45 μ L CHClCOCl (98%) to a 20 mL flask, and dilute it to 20 mL.

Organic impurities mixture standard correction: addition of 150, 180, and 200 μ L standard stock to a 20 ml flask respectively, and then diluted to 20 ml. The content of CH₂ClCOCl, CHCl₂COCl, and CCl₃COCl are 0.857×10⁻³ ~ 2.143×10⁻², 1.714×10⁻³ ~ 4.286×10⁻³, 0.952×10⁻⁴ ~ 2.381×10⁻³ mg/g. These solutions were used as calibration.

Determination of organic impurities: The background was scanned while dry nitrogen flow through it at one end and seal the other end with PTFE; then the air tightness was checked by the soap bubble. Injection of 7 mL of organic impurities mixture standard correction and $TiCl_4$ to a vacuum IR liquid cell and the infrared spectrum was determined. The mass fraction of organic impurities in $TiCl_4$ was

calculated by standard addition or spiking method.

2.4. Air tightness test

One end of assembled infrared pool sealed with soft silicone rubber, and dry nitrogen was poured into the pool from the other end, then checked tightness with soap bubbles.

2.5. Overlapping peak resolve

The peak of $CH_{3-n}Cl_nCOCl$ was processed by OMNIC software. Set all the parameters and choose the peak position. Then the peaks were separated and it was convenient to quantify the three organic compounds.

3. Results and discussion

3.1. Stability test of IR liquid cell

An organic impurities standard calibration was added to IR liquid cell. And we recorded the absorbance of CCl_3COCl , $CHCl_2COCl$ and $CH_2ClCOCl$ in vacuum space every day for 7 days, and then calculate the content of each impurity. The determination result showed that no obvious difference on organic impurities content with the increase of storage period, and all RSD is less than 10%. It indicated this infrared liquid cell presented a superior inertia and seal.

3.2. Precision test

Table 1 indicated the content of $CH_2ClCOCl$, $CHCl_2COCl$ and CCl_3COCl in $TiCl_4$, and each sample was tested for five times. When RSD was less than 0.04%, it indicated a good precision of IR liquid cell which can meet the study requirement.

Table 1. The results of the same organic impurities mixed standard sample at five measuring times $(\mu g \cdot g^{-1})$.

Measure time	X (CH ₂ ClCOCl)	X (CHCl ₂ COCl)	X (CCl ₃ COCl)
1	1.7860×10 ⁻²	2.1858×10^{-3}	3.629×10 ⁻³
2	1.7858×10^{-2}	2.1857×10 ⁻³	3.628×10 ⁻³
3	1.7859×10^{-2}	2.1859×10 ⁻³	3.626×10 ⁻³
4	1.7859×10^{-2}	2.1859×10 ⁻³	3.626×10 ⁻³
5	1.7857×10^{-2}	2.1857×10 ⁻³	3.627×10 ⁻³
RSD (%)	0.0064	0.0045	0.0359

3.3. Detection limit by standard addition

Calculate detection limit by the lowest content of $CHCl_2COCl$ in $TICl_4$ (absorbance=0.003). The calculated detection limit of CCl_3COCl , $CHCl_2COCl$ and $CH_2ClCOCl$ were 3.159×10^{-3} , 1.917×10^{-3} and 1.554×10^{-2} mg/g, respectively. The measured value of CCl_3COCl , $CHCl_2COCl$ and $CH_2ClCOCl$ were lower than level standard of the compounds in fine $TiCl_4$ in Russia. So the precision can meet the quality control and analysis requirement of impurities in $TiCl_4$ production.

3.4. Sample determination

The content of CCl₃COCl, CHCl₂COCl and CH₂ClCOCl in TiCl₄ were tested in IR liquid cell with a standard. Figure 2(a) showed peaks of CCl₃COCl, CHCl₂COCl and CH₂ClCOCl were highly overlapped, which is not able to quantify directly. Carbonyl group (C-O) presented at CCl₃COCl, CHCl₂COCl and CH₂ClCOCl with a sensitive peak at around 1800 cm⁻¹. Inductive effect and other synergistic effect to carbonyl group occurred and it varied from the number of chlorine atoms attached to methyl group. The effects made the difference of sensitive absorption peak, however, overlapping presented when the three organic was determined simultaneously. So the peaks were then separated by OMNIC software [13] for quantifying and the figure was showed in figure 2(b). The infrared

spectroscopy characteristic peak of CCl₃COCl, CHCl₂COCl and CH₂ClCOCl are 1801, 1809, and 1820 cm⁻¹ respectively. The sample determination result was illustrated in table 2. These formulas are regression equation of the three impurity compound. The X illustrated the content of the impurity and the mean of Y was measure value of the compound by FT-IR. The correlation coefficient is 0.8995 for CH₂ClCOCl standard curve, which is the only one less than 0.996. The average level of CH₂ClCOCl, CHCl₂COCl and CCl₃COCl in TiCl₄ was 1.7859×10^{-2} , 2.1857×10^{-3} and 3.626×10^{-3} µg/g respectively. Furthermore, the correlation among the chloroacetic acids and their wavelengths of maximum absorbance was also tested. But no obvious correlation (linear or quadratics) can be found through the result.




Figure 2. Infrared spectra of the organic impurities in $TiCl_4((a)$ showed the original results by FT-IR, (b) showed the results processed by OMNIC software).

Compound	Regression equation	Correlation coefficient(R ²)	Content($\mu g \cdot g^{-1}$)
CH ₂ ClCOCl	$y=-1.571x\times10^{-5}+1.00001785946$	0.8995	1.7859×10^{-2}
CHCl ₂ COCl	$y=-4.205x\times10^{-4}+1.00000223545$	0.9999	2.1857×10^{-3}
CCl ₃ COCl	$y=-1.273x\times10^{-5}+1.0000036264$	0.9962	3.626×10 ⁻³

Table 2. The concentration of the three organic impurities in TiCl₄.

4. Conclusions

A new device with a window film of ZnSe and a body of PTFE can be used to determine the content of impurities in TiCl₄. The regression equation of CH₂ClCOCl was y=-1.571x×10⁻⁵+1.00001785946, R²=0.8995; the regression equation of CHCl₂COCl was y=-4.205x×10⁻⁴+1.00000223545, R²=0.9999; and the regression equation of CCl₃COCl was y=-1.273x×10⁻⁵+1.0000036264, R²=0.9963. The new infrared measuring device showed various advantages, including safe, simple, easily removable, washable and reusable. The process of peak by OMNIC software made the research convenient due to the measurement of different compound simultaneously. This new method prevented the influence of added solution on the matrix, which is applicable for multi-component analysis.

Acknowledgments

The research of this paper was financially supported by the Youth Fund in Academy of Guizhou Science. The test was done in Guizhou Academy of testing and analysis with FT-IR. In developing the ideas presented here, I have received helpful input from Tan H and He J L. I have to thank my co-author, Wang Y J.

References

- [1] Deng G Z and Wang X D 2007 Promote technological progress is key to maintaining sustainable development of titanium sponge production *Development of Titanium Production* 24 9-14
- [2] Ouyang Q S, Zhao Z W and Zhu Y H 2004 Advances in technology and the new titanium sponge production Rare Metals and Cemented Carbides **32** 47-52
- [3] Mo W, Deng G Z and Luo F C 1998 *Titanium Alloys* 1st edn. (Beijing Metallurgical Industry Press, No.39, Zhuyuan North Street, Dongcheng District, Beijing)
- [4] Tan H, He J and Song G L 2012 Determination of trace organic impurities in refined TiCl₄ by infrared spectroscopy *Journal of Instrumental Analysis* **31** 1325-9
- [5] Rolf B and Johnson J 1993 Infrared spectroscopy applied to the determination of impurities in titanium tetrachloride *Technology of Zunyi Titanium* **Sup** 10-5
- [6] Liu W J 1994 Qualitative discussion about infrared spectra of impurities in titanium tetrachloride *Vanadium and Titanium* **21** 26-30
- [7] China Nonferrous Metals Committe (CNSMQ) 2007 Fine titanium tetrachloride quality standards (CN. Q/WBH001-2007)
- [8] Chen W S and Lan W 1996 Determination of fine titanium tetrachloride titanium chloride and oxychloride infrared semi-quantitative method *Guizhou Science* **14** 56-9
- [9] Chen W S and Lan W 2001 Quantitative determination for the trichloro-acetyl chloride in titanium tetrachloride by means of standard addition extrapolation *Journal of Guizhou Normal University (Natural Sciences)* **19** 56-8
- [10] Song Y P, Tan H, Sun J P and He J L 2008 Infrared quantitative analysis of titanium tetrachloride trichloroacetyl chloride *Chinese Determination Technology* **34** 95-7
- [11] Yang L, Huang H and Cheng X L 1999 New calculation method of standard addition method in infrared spectroscopy quantitative analysis *Analytical Laboratory* **18** 79-82
- [12] Griffiths P R and Haseth J A D E 2007 Fourier Transform Infrared Spectroscopy 2nd edn. (Hoboken: John Wiley & Sons Ltd)
- [13] Cai X L 2000 Identification of infrared difference spectroscopy for mixing hypnotics Spectrosc and Spectral Analysis 20 329-32